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Function and indication

All IST® appliances from Prof. Hinz have **two functions**:

1. Fix the lower jaw during sleep in a protruded position and
2. at the same time, reposition the tongue, tongue base and adjacent soft tissue holding all in forward position.

This helps to:

- achieve an enlargement of the rear pharyngeal space
- prevent the rear airway from being occluded or obstructed by a collapsing tongue.

The IST® appliance is designed to:

- reduce or avoid primary – harmless – snoring, and
- prevent temporary respiratory arrest in cases of a mild obstructive sleep apnea syndrome.

Advantages of the IST® appliance:

- thin thermoformed splints do not restrict tongue movement
- unlimited mouth opening movements
- undisturbed mouth breathing
- lateral movements are possible
- good retention without wire elements
- individual adjustment of the lower jaw forward displacement by the dentist 1-2 mm

1. Occlusion of rear airway by the collapsing tongue.


2. Mechanical forward displacement of the lower jaw by a IST® appliance including the tongue


Preparations

1. In order to obtain the master models, these should be duplicated or the existing one impression should be cast 2x.
2. The working models shall be articulated with the aid of the IST®bite fork in the therapeutic forward position. With the help of the parallelometer, the prosthetic equator of each tooth is determined and marked.
3. Undertrained sites such as pontics, switching gaps or interdental spaces are to be blocked out with wax.
4. Before all thermoforming operations, when embedding the models in the granules, care should be taken that a distance of about 3 to 4 mm from the sulcus to the granulate is maintained.
5. The insulating film is deep drawn according to the data. In advance of the subsequent deep drawing, it is advisable to cut the insulating film interdentally by means of a scalpel and trim the film around the model.
6. A 2 mm thick film is now pulled over the prepared model e.g. either Duran or LAMlone.



The IST[®] classic appliance

- 1 Deep draw splints for the upper and lower jaw
- 2 Telescopes
- 3 Incorporated fixing elements
- 4 Protrusion screw for infinitely variable adjustment
- 5 Lateral bite planes in the upper jaw and lower jaw
- 6 Screw for fixing telescopes

Placement of the telescopes

1. Next, the orientation and marking of the installation areas takes place for the telescopes. In order to enable a subsequent adjustment and, as the case may be, distally as well, the telescopes should be pre-activated by 2 mm before the installation.
2. The attachments in the upper jaw are interocclusively and vestibularly in the vicinity of the molars. In the lower jaw, the fixing elements are incorporated interocclusively and vestibularly in the vicinity of the first and/or second premolars. The fixing elements should be at right angles to the middle of the alveolar ridge.
3. The surfaces of the splints and the retainers should be blasted before fixing (in order to prevent crevice corrosion).

Manufacturing of the bite blocks

The upper jaw bite blocks are made of acrylate in the articulator. They are exactly opposite the upper jaw fastening parts in the vicinity of the eye tooth. Check the fit of the rods. The lower jaw bites blocks are located in the vicinity of the lower jaw fastening parts and accordingly, they have to be adapted to the upper jaw bite blocks.

Completion of the IST® classic appliance

The connection of the telescopes to the fastening parts takes place using screws (these can be fixed with screw locking adhesive). The HF 15 is used when the HF 9 is insufficient.

Warning note: The fastening parts should not be damaged when finishing and polishing outside

The IST® classic appliance should be cleaned daily with a toothbrush under cold water. As the case may be, use our **cleaning tabs** (Art. no. 95 010) and store in a dry place.

The IST[®] classic^{NEW} appliance

- 1 Deep draw splints
- 2 IST[®]classic^{NEW} guiding bridges
- 3 Threaded rods (individually and continuously adjustable)
- 4 Lateral bite block in the upper jaw and lower jaw

WARNING: New fitting position!

Adjustment of the bridge supports

1. Put the working models with film drawn over, back into the articulator.
2. In the upper jaw, the fastening parts are in the vicinity of the eye tooth. In the lower jaw, they are in the vicinity of the molars.
3. Fix the fastening parts with Sil-Tech putty and position with the help of the rods. Pay attention to the parallelism of the fastening parts.
4. It may be necessary to adapt the retentions and/or occlusion to the tooth shape.
5. They are fixed with acrylate.

In order to continue protruding after completing the device, the rod threads can be shortened. This is the task of the technician.

WARNING: New fitting position!

Manufacturing of the bite blocks

The lower jaw bite blocks are made from acrylate in the articulator. Here, bite blocks are only required in the eye tooth/premolar vicinity. Check the fit of the rods. Upper jaw bite blocks lie on the upper jaw fastening part and have to be adjusted to the lower jaw.

Protrude, but how?

A change to the protrusion should only be done by the dentist providing the treatment. In order to change the setting on the appliance later, please note the following: An Allen screw must be removed on each side so that the rod can be shortened/extended. Then replace the Allen screw. **The rod must be the same length on both sides!**

Important: **Protrusion** = Shortening the threaded rod (twisting in)
Retrusion = Extending the threaded rod (twisting out)

The IST®classic^{NEW} appliance should be washed on a daily basis with a toothbrush under cold water. If necessary, our **cleaning tabs** (Art. no. 95 010) should be used. Store in a dry place.

Notes on processing the splints:

Before lifting, mark the splint extension and prepare the splints as usual.


The finished lower jaw splint.

The ISTplus® appliance

- 1 Adjustable connecting bridge
- 2 Adjustment screw
- 3 Bridge screw with Allen screw
- 4 Deep drawing splints for upper and lower jaw
- 5 Plastic anchoring part

- 1 Adjustable connecting bridge
- 2 Adjustment screw
- 3 Bridge screw with Allen screw
- 4 Space saver
- 5 Bridge nut
- 6 Plastic anchoring part

Placement of the connecting bridges

1. The adjustable connecting bars are positioned in such a way, that they point mesially in the upper jaw and distally in the lower jaw. The Allen key of the screw should point buccally. The complete plastic anchoring parts and the basis of the splint are roughened, wherein the sides of the anchoring parts are roughened from the outside.
2. It should be noted that the coloured space holders are in the plastic anchoring parts and that they have contact with deep drawing film. It is only this way that it is guaranteed that the connecting bridges are mounted to the maximum extent and at the same time, no plastic gets to the bridge screw. The plastic anchoring parts should run as horizontally as possible to the middle of the alveolar ridge. Now the plastic anchoring parts are fixed to the base with some Sil-Tech putty and in addition, they are fixed in the lower jaw in the thinly prepared plastic. By analogy, the plastic anchoring parts are fastened in the upper jaw.
3. Loosen the bridge screws using the Allen key and remove the connecting bridge.
4. In the lower jaw, small plastic bites blocks are attached. Its position and orientation depends on the OK plastic anchors.
5. The splints are detached from the model, separated below the marked equator and finished. Then remove the coloured space holders from the plastic anchoring parts and the insulating foil from the splints. A light polishing or coating takes place with a light-polymerising paint (PLAQUIT / Dreve).

Completion of the ISTplus® appliance

The connecting bridges are finally mounted with the bridge screws in the plastic anchoring parts.

Important: The mounting of the connecting bridges must be done at the end so that the laser marking on the connecting bridge points mesially on the right hand side and distally on the left hand side.

In order to change the protrusion, please note the following: A **shortening** of the connecting bridge reinforces the protrusion, but the adjustment screw is turned from the upper jaw to the lower jaw (downwards). In order to reduce the protrusion, the connecting bridge has to be extended by turning the adjusting screw from the lower jaw to the upper jaw (upwards). One revolution (360 °) provides a change of 0.7 mm.

The ISTplus® appliance should be washed daily with a toothbrush under cold water. If necessary, use our **cleaning tabs** (Art. no. 95 010) and store in a dry place.

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